

配对 抗体

Matched
Antibody Pair

高特异性

高信噪比

包含常见细胞因子、肿瘤
或其他疾病关键标志物等

抗体特异性经过
阳性细胞裂解液验证

可用于夹心法 ELISA
和其他免疫检测



成都正能生物技术有限责任公司



PlatinumAb
精品抗体



RecRAbs®
重组兔单抗



Hist-Plus®
病理级抗体



Ab Pairs
抗体对




KnockOut
KO验证抗体

微信公众号




ELISA 检测原理示意图



配对抗体 (Matched Antibody Pair) 指可以同时结合在一个抗原分子上识别不同表位的两个抗体, 一般应用于双抗夹心法 ELISA, 流式液相多重检测 (Cytometric Bead Array, CBA), 也可能适用于侧向层析 (LFA) 或其他相关的免疫检测。

正能提供的抗体对为小鼠单克隆抗体或重组兔单克隆抗体配对, 均经过 ELISA 检测验证。抗体对产品线覆盖细胞因子、免疫球蛋白、肿瘤标志物或其他疾病相关蛋白等。我们还可提供匹配的抗体对提供蛋白质标准品作为对照或生成标准曲线。

部分产品的 ELISA 验证图



产品信息说明

货号前缀

MS: 抗体来源小鼠
RB: 抗体来源为兔
RT: 抗体来源为大鼠

货号后缀

UC: 未标记的捕获抗体
UD: 未标记的检测抗体
HD: 带 HRP 标记的检测抗体

产品名称后缀

capture: 捕获抗体
detector: 检测抗体

靶标	货号	产品名称	反应种属	应用
AFP	MS20038UC	alpha 1 Fetoprotein Mouse mAb, capture	Human	ELISA
	MS20038UD	alpha 1 Fetoprotein Mouse mAb, unconjugated, detector	Human	ELISA
Alpha-synuclein	RB20116UC	Alpha-synuclein Rabbit mAb, capture	Human	ELISA
	RB20116UD	Alpha-synuclein Rabbit mAb, unconjugated, detector	Human	ELISA
AMBP	MS20010UC	AMBP Mouse mAb, capture	Human	ELISA
	MS20010UD	AMBP Mouse mAb, unconjugated, detector	Human	ELISA
BACE1	RB20119UC	BACE1 Rabbit mAb, capture	Human	ELISA
	RB20119UD	BACE1 Rabbit mAb, unconjugated, detector	Human	ELISA
CCL18	RB20063UC	CCL18 Rabbit mAb, capture	Human	ELISA
	RB20063UD	CCL18 Rabbit mAb, unconjugated, detector	Human	ELISA
CD147	RB20124UC	CD147 Rabbit mAb, capture	Human	ELISA
	RB20124UD	CD147 Rabbit mAb, unconjugated, detector	Human	ELISA
CD40	RB20126UC	CD40 Rabbit mAb, capture	Human	ELISA
	RB20126UD	CD40 Rabbit mAb, unconjugated, detector	Human	ELISA
CD68	RB20128UC	CD68 Rabbit mAb, capture	Human	ELISA
	RB20128UD	CD68 Rabbit mAb, unconjugated, detector	Human	ELISA
CEA	MS20030UC	CEA Mouse mAb, capture	Human	ELISA
	MS20030UD	CEA Mouse mAb, unconjugated, detector	Human	ELISA
CGB3	MS20011UC	CGB3 Mouse mAb, capture	Human	ELISA
	MS20011UD	CGB3 Mouse mAb, unconjugated, detector	Human	ELISA
COMP	RB20129UC	COMP/Cartilage oligomeric matrix protein Rabbit mAb, capture	Human	ELISA
	RB20129UD	COMP/Cartilage oligomeric matrix protein Rabbit mAb, unconjugated, detector	Human	ELISA
CTLA4	MS20012UC	CTLA4 Mouse mAb, capture	Human	ELISA
	MS20012UD	CTLA4 Mouse mAb, unconjugated, detector	Human	ELISA
CXCL9	RB20066UC	CXCL9 Rabbit mAb, capture	Human	ELISA
	RB20066UD	CXCL9 Rabbit mAb, unconjugated, detector	Human	ELISA
Cystatin C	MS20046UC	Cystatin C Mouse mAb, capture	Human	ELISA
	MS20046UD	Cystatin C Mouse mAb, unconjugated, detector	Human	ELISA
Cytokeratin 19	RB20062UC	Cytokeratin 19 Rabbit mAb, capture	Human	ELISA
	RB20062UD	Cytokeratin 19 Rabbit mAb, unconjugated, detector	Human	ELISA
CHI3L1	RB20201UC	YKL-40/CHI3L1 Rabbit mAb, capture	Human	ELISA
	RB20201UD	YKL-40/CHI3L1 Rabbit mAb, unconjugated, detector	Human	ELISA
DCX	RB20130UC	Doublecortin Rabbit mAb, capture	Human	ELISA
	RB20130UD	Doublecortin Rabbit mAb, unconjugated, detector	Human	ELISA
E Cadherin	RB20133UC	E Cadherin Rabbit mAb, capture	Human	ELISA
	RB20133UD	E Cadherin Rabbit mAb, unconjugated, detector	Human	ELISA
EGF	MS20035UC	EGF Mouse mAb, capture	Human	ELISA
	MS20035UD	EGF Mouse mAb, unconjugated, detector	Human	ELISA
	MS20048UC	EGF Mouse mAb, capture	Human	ELISA
	MS20048UD	EGF Mouse mAb, unconjugated, detector	Human	ELISA
Endostatin	MS20013UC	Endostatin Mouse mAb, capture	Human	ELISA
	MS20013UD	Endostatin Mouse mAb, unconjugated, detector	Human	ELISA
ENO1	MS20031UC	ENO1 Mouse mAb, capture	Human	ELISA
	MS20031UD	ENO1 Mouse mAb, unconjugated, detector	Human	ELISA
EpCAM	RB20115UC	EpCAM Rabbit mAb, capture	Human	ELISA
	RB20115UD	EpCAM Rabbit mAb, unconjugated, detector	Human	ELISA
FABP	MS20052UC	FABP1 Mouse mAb, capture	Human	ELISA
	MS20052UD	FABP1 Mouse mAb, unconjugated, detector	Human	ELISA
	MS20047UC	FABP3 Mouse mAb, capture	Human	ELISA
	MS20047UD	FABP3 Mouse mAb, unconjugated, detector	Human	ELISA
Factor B	MS20014UC	Factor B Mouse mAb, capture	Human	ELISA
	MS20014UD	Factor B Mouse mAb, unconjugated, detector	Human	ELISA
FGF2	MS20015UC	FGF2 Mouse mAb, capture	Human	ELISA
	MS20015UD	FGF2 Mouse mAb, unconjugated, detector	Human	ELISA
	RB20015UC	FGF2 Rabbit mAb, capture	Human	ELISA
	RB20015UD	FGF2 Rabbit mAb, unconjugated, detector	Human	ELISA
FGA	RT20037UC	Fibrinogen alpha Chain Rat mAb, capture	Human	ELISA
	RT20037UD	Fibrinogen alpha Chain Rat mAb, unconjugated, detector	Human	ELISA
FOXA2	RB20134UC	FOXA2 Rabbit mAb, capture	Human	ELISA
	RB20134UD	FOXA2 Rabbit mAb, unconjugated, detector	Human	ELISA
GAD	RB20135UC	GAD65 Rabbit mAb, capture	Human	ELISA
	RB20135UD	GAD65 Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20136UC	GAD67 Rabbit mAb, capture	Human	ELISA
	RB20136UD	GAD67 Rabbit mAb, unconjugated, detector	Human	ELISA
GAP43	RB20138UC	GAP43 Rabbit mAb, capture	Human	ELISA
	RB20138UD	GAP43 Rabbit mAb, unconjugated, detector	Human	ELISA

靶标	货号	产品名称	反应种属	应用
G-CSF	MS20016UC	G-CSF Mouse mAb, capture	Human	ELISA
	MS20016UD	G-CSF Mouse mAb, unconjugated, detector	Human	ELISA
	RB20016UC	G-CSF Rabbit mAb, capture	Human	ELISA
	RB20016UD	G-CSF Rabbit mAb, unconjugated, detector	Human	ELISA
GFAP	RB20139UC	GFAP Rabbit mAb, capture	Human	ELISA
	RB20139UD	GFAP Rabbit mAb, unconjugated, detector	Human	ELISA
GM-CSF	RB20017UC	GM-CSF Rabbit mAb, capture	Human	ELISA
	RB20017UD	GM-CSF Rabbit mAb, unconjugated, detector	Human	ELISA
	RT20017UC	GM-CSF Rat mAb, capture	Human	ELISA
	RT20017UD	GM-CSF Rat mAb, unconjugated, detector	Human	ELISA
CXCL1	RB20064UC	GRO alpha Rabbit mAb, capture	Human	ELISA
	RB20064UD	GRO alpha Rabbit mAb, unconjugated, detector	Human	ELISA
GHI	RT20034UC	Growth hormone 1 Rat mAb, capture	Human	ELISA
	RT20034UD	Growth hormone 1 Rat mAb, unconjugated, detector	Human	ELISA
HGF	RB20067UC	Hepatocyte growth factor Rabbit mAb, capture	Human	ELISA
	RB20067UD	Hepatocyte growth factor Rabbit mAb, unconjugated, detector	Human	ELISA
IFN alpha	RB20068UC	IFN alpha1 Rabbit mAb, capture	Human	ELISA
	RB20068UD	IFN alpha1 Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20089UC	IFN alpha2 Rabbit mAb, capture	Human	ELISA
	RB20089UD	IFN alpha2 Rabbit mAb, unconjugated, detector	Human	ELISA
IFN gamma	MS20026UC	IFN gamma Mouse mAb, capture	Human	ELISA
	MS20026UD	IFN gamma Mouse mAb, unconjugated, detector	Human	ELISA
IgA	MS20018UC	IgA Mouse mAb, capture	Human	ELISA
	MS20018UD	IgA Mouse mAb, unconjugated, detector	Human	ELISA
IgE	MS20019UC	IgE Mouse mAb, capture	Human	ELISA
	MS20019UD	IgE Mouse mAb, unconjugated, detector	Human	ELISA
IgG	MS20042UC	IgG Mouse mAb, capture	Human	ELISA
	MS20042UD	IgG Mouse mAb, unconjugated, detector	Human	ELISA
	RT20057UC	IgG1 Rat mAb, capture	Human	ELISA
	MS20057UD	IgG Mouse mAb, unconjugated, detector	Human	ELISA
	MS20058UC	IgG2 Mouse mAb, capture	Human	ELISA
	MS20058UD	IgG2 Mouse mAb, unconjugated, detector	Human	ELISA
	MS20059UC	IgG3 Mouse mAb, capture	Human	ELISA
	MS20059UD	IgG Mouse mAb, unconjugated, detector	Human	ELISA
	MS20060UC	IgG4 Mouse mAb, capture	Human	ELISA
	MS20060UD	IgG4 Mouse mAb, unconjugated, detector	Human	ELISA
IgM	MS20043UC	IgM Mouse mAb, capture	Human	ELISA
	MS20043UD	IgM Mouse mAb, unconjugated, detector	Human	ELISA
IL-1 alpha	RB20069UC	IL-1 alpha Rabbit mAb, capture	Human	ELISA
	RB20069UD	IL-1 alpha Rabbit mAb, unconjugated, detector	Human	ELISA
IL-1 beta	MS20040UC	IL-1 beta Mouse mAb, capture	Human	ELISA
	MS20040UD	IL-1 beta Mouse mAb, unconjugated, detector	Human	ELISA
	RB20040UC	IL-1 beta Rabbit mAb, capture	Human	ELISA
	RB20040UD	IL-1 beta Rabbit mAb, unconjugated, detector	Human	ELISA
IL-1RA	RB20079UC	IL-1RA Rabbit mAb, capture	Human	ELISA
	RB20079UD	IL-1RA Rabbit mAb, unconjugated, detector	Human	ELISA
IL-2	MS20027UC	IL-2 Mouse mAb, capture	Human	ELISA
	MS20027UD	IL-2 Mouse mAb, unconjugated, detector	Human	ELISA
IL-2RA	RB20081UC	IL-2RA Rabbit mAb, capture	Human	ELISA
	RB20081UD	IL-2RA Rabbit mAb, unconjugated, detector	Human	ELISA
IL-3	RB20082UC	IL-3 Rabbit mAb, capture	Human	ELISA
	RB20082UD	IL-3 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-4	RB20083UC	IL-4 Rabbit mAb, capture	Human	ELISA
	RB20083UD	IL-4 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-5	RB20084UC	IL-5 Rabbit mAb, capture	Human	ELISA
	RB20084UD	IL-5 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-6	MS20021UC	IL-6 Mouse mAb, capture	Human	ELISA
	MS20021UD	IL-6 Mouse mAb, unconjugated, detector	Human	ELISA
	RB20021UC	IL-6 Rabbit mAb, capture	Human	ELISA
	RB20021UD	IL-6 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-7	RB20085UC	IL-7 Rabbit mAb, capture	Human	ELISA
	RB20085UD	IL-7 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-8	MS20039UC	KC(IL-8) Mouse mAb, capture	Human	ELISA
	MS20039UD	KC(IL-8) Mouse mAb, unconjugated, detector	Human	ELISA
	RB20061UC	IL-8 Rabbit mAb, capture	Human	ELISA
	RB20061UD	IL-8 Rabbit mAb, unconjugated, detector	Human	ELISA

靶标	货号	产品名称	反应种属	应用
IL-9	RB20088UC	IL-9 Rabbit mAb, capture	Human	ELISA
	RB20088UD	IL-9 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-10	RB20070UC	IL-10 Rabbit mAb, capture	Human	ELISA
	RB20070UD	IL-10 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-11	MS20020UC	IL-11 Mouse mAb, capture	Human	ELISA
	MS20020UD	IL-11 Mouse mAb, unconjugated, detector	Human	ELISA
IL-12	RB20071UC	IL-12 p70/IL-12A Rabbit mAb, capture	Human	ELISA
	RB20071UD	IL-12 p70/IL-12A Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20072UC	IL-12B Rabbit mAb, capture	Human	ELISA
	RB20072UD	IL-12B Rabbit mAb, unconjugated, detector	Human	ELISA
IL-13	RB20073UC	IL-13 Rabbit mAb, capture	Human	ELISA
	RB20073UD	IL-13 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-15	RB20074UC	IL-15 Rabbit mAb, capture	Human	ELISA
	RB20074UD	IL-15 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-16	RB20075UC	IL-16 Rabbit mAb, capture	Human	ELISA
	RB20075UD	IL-16 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-17	RB20076UC	IL-17A Rabbit mAb, capture	Human	ELISA
	RB20076UD	IL-17A Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20077UC	IL-17F Rabbit mAb, capture	Human	ELISA
	RB20077UD	IL-17F Rabbit mAb, unconjugated, detector	Human	ELISA
IL-18	RB20078UC	IL-18 Rabbit mAb, capture	Human	ELISA
	RB20078UD	IL-18 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-20	RB20208UC	IL-20 Rabbit mAb, capture	Human	ELISA
	RB20208UD	IL-20 Rabbit mAb, unconjugated, detector	Human	ELISA
IL-22	RB20080UC	IL-22 Rabbit mAb, capture	Human	ELISA
	RB20080UD	IL-22 Rabbit mAb, unconjugated, detector	Human	ELISA
Insulin	RB20209UC	Insulin Rabbit mAb, capture	Human	ELISA
	RB20209UD	Insulin Rabbit mAb, unconjugated, detector	Human	ELISA
Kappa Light Chain	MS20041UC	Kappa Light Chain Mouse mAb, capture	Human	ELISA
	MS20041UD	Kappa Light Chain Mouse mAb, unconjugated, detector	Human	ELISA
KITLG	RT20022UC	Kit ligand Rat mAb, capture	Human	ELISA
	RT20022UD	Kit ligand Rat mAb, unconjugated, detector	Human	ELISA
Lactoferrin	MS20023UC	Lactoferrin Mouse mAb, capture	Human	ELISA
	MS20023UD	Lactoferrin Mouse mAb, unconjugated, detector	Human	ELISA
LIF	RB20090UC	LIF Rabbit mAb, capture	Human	ELISA
	RB20090UD	LIF Rabbit mAb, unconjugated, detector	Human	ELISA
MCP-1	MS20044UC	MCP-1 Mouse mAb, capture	Human	ELISA
	MS20044UD	MCP-1 Mouse mAb, unconjugated, detector	Human	ELISA
	RB20044UC	MCP-1 Rabbit mAb, capture	Human	ELISA
	RB20044UD	MCP-1 Rabbit mAb, unconjugated, detector	Human	ELISA
M-CSF	RB20091UC	M-CSF Rabbit mAb, capture	Human	ELISA
	RB20091UD	M-CSF Rabbit mAb, unconjugated, detector	Human	ELISA
MIF	RB20092UC	MIF Rabbit mAb, capture	Human	ELISA
	RB20092UD	MIF Rabbit mAb, unconjugated, detector	Human	ELISA
MIP-1	RB20093UC	MIP-1 alpha Rabbit mAb, capture	Human	ELISA
	RB20093UD	MIP-1 alpha Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20113UC	MIP-1 beta Rabbit mAb, capture	Human	ELISA
	RB20113UD	MIP-1 beta Rabbit mAb, unconjugated, detector	Human	ELISA
MMP1	RB20142UC	MMP1 Rabbit mAb, capture	Human	ELISA
	RB20142UD	MMP1 Rabbit mAb, unconjugated, detector	Human	ELISA
MMP2	MS20033UC	MMP2 Mouse mAb, capture	Human	ELISA
	MS20033UD	MMP2 Mouse mAb, unconjugated, detector	Human	ELISA
	RB20148UC	MMP2 Rabbit mAb, capture	Human	ELISA
	RB20148UD	MMP2 Rabbit mAb, unconjugated, detector	Human	ELISA
MMP3	RB20151UC	MMP3 Rabbit mAb, capture	Human	ELISA
	RB20151UD	MMP3 Rabbit mAb, unconjugated, detector	Human	ELISA
MMP8	RB20153UC	MMP8 Rabbit mAb, capture	Human	ELISA
	RB20153UD	MMP8 Rabbit mAb, unconjugated, detector	Human	ELISA
MMP9	RB20087UC	MMP9 Rabbit mAb, capture	Human	ELISA
	RB20087UD	MMP9 Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20155UC	MMP9 Rabbit mAb, capture	Human	ELISA
	RB20155UD	MMP9 Rabbit mAb, unconjugated, detector	Human	ELISA
MMP12	RB20145UC	MMP12 Rabbit mAb, capture	Human	ELISA
	RB20145UD	MMP12 Rabbit mAb, unconjugated, detector	Human	ELISA
MMP13	RB20146UC	MMP13 Rabbit mAb, capture	Human	ELISA
	RB20146UD	MMP13 Rabbit mAb, unconjugated, detector	Human	ELISA

靶标	货号	产品名称	反应种属	应用
NPPB	MS20053UC	Natriuretic Peptides B Mouse mAb, capture	Human	ELISA
	MS20053UD	Natriuretic Peptides B Mouse mAb, unconjugated, detector	Human	ELISA
Neurofilament	RB20160UC	Neurofilament Rabbit mAb, capture	Human	ELISA
	RB20160UD	Neurofilament Rabbit mAb, unconjugated, detector	Human	ELISA
Neurogranin	RB20161UC	Neurogranin Rabbit mAb, capture	Human	ELISA
	RB20161UD	Neurogranin Rabbit mAb, unconjugated, detector	Human	ELISA
NG2	RB20163UC	NG2 Rabbit mAb, capture	Human	ELISA
	RB20163UD	NG2 Rabbit mAb, unconjugated, detector	Human	ELISA
NPTX2	RB20165UC	NPTX2 Rabbit mAb, capture	Human	ELISA
	RB20165UD	NPTX2 Rabbit mAb, unconjugated, detector	Human	ELISA
NSE	RB20166UC	NSE Rabbit mAb, capture	Human	ELISA
	RB20166UD	NSE Rabbit mAb, unconjugated, detector	Human	ELISA
P75 Neurogranin	RB20168UC	P75 Neurogranin Rabbit mAb, capture	Human	ELISA
	RB20168UD	P75 Neurogranin Rabbit mAb, unconjugated, detector	Human	ELISA
PDGFR alpha	RB20170UC	PDGFR alpha Rabbit mAb, capture	Human	ELISA
	RB20170UD	PDGFR alpha Rabbit mAb, unconjugated, detector	Human	ELISA
PDGFR beta	RB20171UC	PDGFR beta Rabbit mAb, capture	Human	ELISA
	RB20171UD	PDGFR beta Rabbit mAb, unconjugated, detector	Human	ELISA
PGP9.5	RB20173UC	PGP9.5 Rabbit mAb, capture	Human	ELISA
	RB20173UD	PGP9.5 Rabbit mAb, unconjugated, detector	Human	ELISA
PSA	MS20032UC	PSA Mouse mAb, capture	Human	ELISA
	MS20032UD	PSA Mouse mAb, unconjugated, detector	Human	ELISA
RANTES	RB20094UC	RANTES Rabbit mAb, capture	Human	ELISA
	RB20094UD	RANTES Rabbit mAb, unconjugated, detector	Human	ELISA
RBP4	MS20050UC	RBP4 Mouse mAb, capture	Human	ELISA
	MS20050UD	RBP4 Mouse mAb, unconjugated, detector	Human	ELISA
SNCG	RB20137UC	gamma Synuclein/SNCG Rabbit mAb, capture	Human	ELISA
	RB20137UD	gamma Synuclein/SNCG Rabbit mAb, unconjugated, detector	Human	ELISA
S100A1	RB20178UC	S100A1 Rabbit mAb, capture	Human	ELISA
	RB20178UD	S100A1 Rabbit mAb, unconjugated, detector	Human	ELISA
S100A4	RB20180UC	S100A4 Rabbit mAb, capture	Human	ELISA
	RB20180UD	S100A4 Rabbit mAb, unconjugated, detector	Human	ELISA
S100A9	RB20182UC	S100A9 Rabbit mAb, capture	Human	ELISA
	RB20182UD	S100A9 Rabbit mAb, unconjugated, detector	Human	ELISA
S100A12	RB20179UC	S100A12/CGRP Rabbit mAb, capture	Human	ELISA
	RB20179UD	S100A12/CGRP Rabbit mAb, unconjugated, detector	Human	ELISA
S100B	MS20036UC	S100B Mouse mAb, capture	Human	ELISA
	MS20036UD	S100B Mouse mAb, unconjugated, detector	Human	ELISA
	RB20177UC	S100B Rabbit mAb, capture	Human	ELISA
	RB20177UD	S100B Rabbit mAb, unconjugated, detector	Human	ELISA
SOX2	RB20183UC	SOX2 Rabbit mAb, capture	Human	ELISA
	RB20183UD	SOX2 Rabbit mAb, unconjugated, detector	Human	ELISA
SOX9	RB20184UC	SOX9 Rabbit mAb, capture	Human	ELISA
	RB20184UD	SOX9 Rabbit mAb, unconjugated, detector	Human	ELISA
STMN1	RB20185UC	Stathmin 1 Rabbit mAb, capture	Human	ELISA
	RB20185UD	Stathmin 1 Rabbit mAb, unconjugated, detector	Human	ELISA
TUBB	RB20123UC	beta Tubulin Rabbit mAb, capture	Human	ELISA
	RB20123UD	beta Tubulin Rabbit mAb, unconjugated, detector	Human	ELISA
TUBB3	RB20122UC	beta III Tubulin Rabbit mAb, capture	Human	ELISA
	RB20122UD	beta III Tubulin Rabbit mAb, unconjugated, detector	Human	ELISA
Tau	RB20186UC	Tau Rabbit mAb, capture	Human	ELISA
	RB20186UD	Tau Rabbit mAb, unconjugated, detector	Human	ELISA
Phospho-Tau	MS20054UC	Tau Mouse mAb, capture	Human	ELISA
	MS20054UD	Phospho-Tau (Thr181) Mouse mAb, unconjugated, detector	Human	ELISA
	MS20055UC	Tau Mouse mAb, capture	Human	ELISA
	MS20055UD	Phospho-Tau (Thr217) Mouse mAb, unconjugated, detector	Human	ELISA
	MS20056UC	Tau Mouse mAb, capture	Human	ELISA
	MS20056UD	Phospho-Tau (Thr231) Mouse mAb, unconjugated, detector	Human	ELISA
TDP43	RB20188UC	TDP43 Rabbit mAb, capture	Human	ELISA
	RB20188UD	TDP43 Rabbit mAb, unconjugated, detector	Human	ELISA
THBD	RB20065UC	Thrombomodulin Rabbit mAb, capture	Human	ELISA
	RB20065UD	Thrombomodulin Rabbit mAb, unconjugated, detector	Human	ELISA
TNF alpha	MS20028UC	TNF alpha Mouse mAb, capture	Human	ELISA
	MS20028UD	TNF alpha Mouse mAb, unconjugated, detector	Human	ELISA
	RB20028UC	TNF alpha Rabbit mAb, capture	Human	ELISA
	RB20028UD	TNF alpha Rabbit mAb, unconjugated, detector	Human	ELISA

靶标	货号	产品名称	反应种属	应用
TNF beta	RB20095UC	TNF beta Rabbit mAb, capture	Human	ELISA
	RB20095UD	TNF beta Rabbit mAb, unconjugated, detector	Human	ELISA
Transferrin	MS20051UC	Transferrin Mouse mAb, capture	Human	ELISA
	MS20051UD	Transferrin Mouse mAb, unconjugated, detector	Human	ELISA
TREMI	RB20189UC	TREMI Rabbit mAb, capture	Human	ELISA
	RB20189UD	TREMI Rabbit mAb, unconjugated, detector	Human	ELISA
TREM2	RB20191UC	TREM2 Rabbit mAb, capture	Human	ELISA
	RB20191UD	TREM2 Rabbit mAb, unconjugated, detector	Human	ELISA
VEGF Receptor	RB20192UC	VEGF Receptor 1 Rabbit mAb, capture	Human	ELISA
	RB20192UD	VEGF Receptor 1 Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20194UC	VEGF Receptor 2 Rabbit mAb, capture	Human	ELISA
	RB20194UD	VEGF Receptor 2 Rabbit mAb, unconjugated, detector	Human	ELISA
	RB20197UC	VEGF Receptor 3 Rabbit mAb, capture	Human	ELISA
	RB20197UD	VEGF Receptor 3 Rabbit mAb, unconjugated, detector	Human	ELISA
Vimentin	RB20199UC	Vimentin Rabbit mAb, capture	Human	ELISA
	RB20199UD	Vimentin Rabbit mAb, unconjugated, detector	Human	ELISA
Alpha-synuclein	RB20117UC	Alpha-synuclein Rabbit mAb, capture	Mouse	ELISA
	RB20117UD	Alpha-synuclein Rabbit mAb, unconjugated, detector	Mouse	ELISA
ANG2	RB20118UC	Angiopoietin 2/ANG2 Rabbit mAb, capture	Mouse	ELISA
	RB20118UD	Angiopoietin 2/ANG2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
CD147	RB20125UC	CD147 Rabbit mAb, capture	Mouse	ELISA
	RB20125UD	CD147 Rabbit mAb, unconjugated, detector	Mouse	ELISA
CD40	RB20127UC	CD40 Rabbit mAb, capture	Mouse	ELISA
	RB20127UD	CD40 Rabbit mAb, unconjugated, detector	Mouse	ELISA
CHI3L1	RB20200UC	YKL-40/CHI3L1 Rabbit mAb, capture	Mouse	ELISA
	RB20200UD	YKL-40/CHI3L1 Rabbit mAb, unconjugated, detector	Mouse	ELISA
E Cadherin	RB20131UC	E Cadherin Rabbit mAb, capture	Mouse	ELISA
	RB20131UD	E Cadherin Rabbit mAb, unconjugated, detector	Mouse	ELISA
GFAP	RB20140UC	GFAP Rabbit mAb, capture	Mouse	ELISA
	RB20140UD	GFAP Rabbit mAb, unconjugated, detector	Mouse	ELISA
IFN γ	RB20106UC	IFN gamma Rabbit mAb, capture	Mouse	ELISA
	RB20106UD	IFN gamma Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-1 beta	RB20096UC	IL-1 beta Rabbit mAb, capture	Mouse	ELISA
	RB20096UD	IL-1 beta Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-2	RB20101UC	IL-2 Rabbit mAb, capture	Mouse	ELISA
	RB20101UD	IL-2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-4	RB20102UC	IL-4 Rabbit mAb, capture	Mouse	ELISA
	RB20102UD	IL-4 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-5	RB20103UC	IL-5 Rabbit mAb, capture	Mouse	ELISA
	RB20103UD	IL-5 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-6	RB20104UC	IL-6 Rabbit mAb, capture	Mouse	ELISA
	RB20104UD	IL-6 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-8	RB20105UC	IL-8 Rabbit mAb, capture	Mouse	ELISA
	RB20105UD	IL-8 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-10	RB20097UC	IL-10 Rabbit mAb, capture	Mouse	ELISA
	RB20097UD	IL-10 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-12	RB20098UC	IL-12 p40 Rabbit mAb, capture	Mouse	ELISA
	RB20098UD	IL-12 p40 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-13	RB20099UC	IL-13 Rabbit mAb, capture	Mouse	ELISA
	RB20099UD	IL-13 Rabbit mAb, unconjugated, detector	Mouse	ELISA
IL-17	RB20100UC	IL-17A Rabbit mAb, capture	Mouse	ELISA
	RB20100UD	IL-17A Rabbit mAb, unconjugated, detector	Mouse	ELISA
MIF	RB20141UC	MIF Rabbit mAb, capture	Mouse	ELISA
	RB20141UD	MIF Rabbit mAb, unconjugated, detector	Mouse	ELISA
MMP2	RB20147UC	MMP2 Rabbit mAb, capture	Mouse	ELISA
	RB20147UD	MMP2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
MMP3	RB20150UC	MMP3 Rabbit mAb, capture	Mouse	ELISA
	RB20150UD	MMP3 Rabbit mAb, unconjugated, detector	Mouse	ELISA
MMP8	RB20154UC	MMP8 Rabbit mAb, capture	Mouse	ELISA
	RB20154UD	MMP8 Rabbit mAb, unconjugated, detector	Mouse	ELISA
MMP9	RB20156UC	MMP9 Rabbit mAb, capture	Mouse	ELISA
	RB20156UD	MMP9 Rabbit mAb, unconjugated, detector	Mouse	ELISA
MMP10	RB20143UC	MMP10 Rabbit mAb, capture	Mouse	ELISA
	RB20143UD	MMP10 Rabbit mAb, unconjugated, detector	Mouse	ELISA
MMP12	RB20144UC	MMP12 Rabbit mAb, capture	Mouse	ELISA
	RB20144UD	MMP12 Rabbit mAb, unconjugated, detector	Mouse	ELISA

靶标	货号	产品名称	反应种属	应用
Mog	RB20157UC	Myelin oligodendrocyte glycoprotein Rabbit mAb, capture	Mouse	ELISA
	RB20157UD	Myelin oligodendrocyte glycoprotein Rabbit mAb, unconjugated, detector	Mouse	ELISA
N Cadherin	RB20158UC	N Cadherin Rabbit mAb, capture	Mouse	ELISA
	RB20158UD	N Cadherin Rabbit mAb, unconjugated, detector	Mouse	ELISA
Nestin	RB20159UC	Nestin Rabbit mAb, capture	Mouse	ELISA
	RB20159UD	Nestin Rabbit mAb, unconjugated, detector	Mouse	ELISA
NG2	RB20162UC	NG2 Rabbit mAb, capture	Mouse	ELISA
	RB20162UD	NG2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
NMDAR1	RB20164UC	NMDAR1 Rabbit mAb, capture	Mouse	ELISA
	RB20164UD	NMDAR1 Rabbit mAb, unconjugated, detector	Mouse	ELISA
NSE	RB20167UC	NSE Rabbit mAb, capture	Mouse	ELISA
	RB20167UD	NSE Rabbit mAb, unconjugated, detector	Mouse	ELISA
p75 NGF Receptor	RB20169UC	p75 NGF Receptor Rabbit mAb, capture	Mouse	ELISA
	RB20169UD	p75 NGF Receptor Rabbit mAb, unconjugated, detector	Mouse	ELISA
PDGFR beta	RB20172UC	PDGFR beta Rabbit mAb, capture	Mouse	ELISA
	RB20172UD	PDGFR beta Rabbit mAb, unconjugated, detector	Mouse	ELISA
PGP9.5	RB20174UC	PGP9.5 Rabbit mAb, capture	Mouse	ELISA
	RB20174UD	PGP9.5 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RAGE	RB20176UC	RAGE Rabbit mAb, capture	Mouse	ELISA
	RB20176UD	RAGE Rabbit mAb, unconjugated, detector	Mouse	ELISA
S100A4	RB20181UC	S100A4 Rabbit mAb, capture	Mouse	ELISA
	RB20181UD	S100A4 Rabbit mAb, unconjugated, detector	Mouse	ELISA
Tau	RB20187UC	Tau Rabbit mAb, capture	Mouse	ELISA
	RB20187UD	Tau Rabbit mAb, unconjugated, detector	Mouse	ELISA
TNF alpha	RB20107UC	TNF alpha Rabbit mAb, capture	Mouse	ELISA
	RB20107UD	TNF alpha Rabbit mAb, unconjugated, detector	Mouse	ELISA
TREM2	RB20190UC	TREM2 Rabbit mAb, capture	Mouse	ELISA
	RB20190UD	TREM2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
VEGF	RB20108UC	VEGFA Rabbit mAb, capture	Mouse	ELISA
	RB20108UD	VEGFA Rabbit mAb, unconjugated, detector	Mouse	ELISA
	RB20198UC	VEGFB Rabbit mAb, capture	Mouse	ELISA
	RB20198UD	VEGFB Rabbit mAb, unconjugated, detector	Mouse	ELISA
VEGF Receptor	RB20195UC	VEGF Receptor 2 Rabbit mAb, capture	Mouse	ELISA
	RB20195UD	VEGF Receptor 2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
	RB20196UC	VEGF Receptor 3 Rabbit mAb, capture	Mouse	ELISA
	RB20196UD	VEGF Receptor 3 Rabbit mAb, unconjugated, detector	Mouse	ELISA
E Cadherin	RB20132UC	E Cadherin Rabbit mAb, capture	Rat	ELISA
	RB20132UD	E Cadherin Rabbit mAb, unconjugated, detector	Rat	ELISA
IgM	MS20024UC	IgM Mouse mAb, capture	Rat	ELISA
	MS20024UD	IgM Mouse mAb, unconjugated, detector	Rat	ELISA
IL-1 beta	RB20109UC	IL-1 beta Rabbit mAb, capture	Rat	ELISA
	RB20109UD	IL-1 beta Rabbit mAb, unconjugated, detector	Rat	ELISA
IL-6	RB20111UC	IL-6 Rabbit mAb, capture	Rat	ELISA
	RB20111UD	IL-6 Rabbit mAb, unconjugated, detector	Rat	ELISA
IL-10	RB20110UC	IL-10 Rabbit mAb, capture	Rat	ELISA
	RB20110UD	IL-10 Rabbit mAb, unconjugated, detector	Rat	ELISA
Kappa Light Chain	MS20025UC	Kappa light chain Mouse mAb, capture	Rat	ELISA
	MS20025UD	Kappa light chain Mouse mAb, unconjugated, detector	Rat	ELISA
MMP3	RB20152UD	MMP3 Rabbit mAb, unconjugated, detector	Rat	ELISA
	RB20152UC	MMP3 Rabbit mAb, capture	Rat	ELISA
RAGE	RB20175UC	RAGE Rabbit mAb, capture	Rat	ELISA
	RB20175UD	RAGE Rabbit mAb, unconjugated, detector	Rat	ELISA
TNF alpha	RB20112UC	TNF alpha Rabbit mAb, capture	Rat	ELISA
	RB20112UD	TNF alpha Rabbit mAb, unconjugated, detector	Rat	ELISA
VEGF Receptor	RB20193UC	VEGF Receptor 1 Rabbit mAb, capture	Rat	ELISA
	RB20193UD	VEGF Receptor 1 Rabbit mAb, unconjugated, detector	Rat	ELISA
Lipocalin 2	MS20049UC	Lipocalin 2 Mouse mAb, capture	Human,Rat	ELISA
	MS20049UD	Lipocalin 2 Mouse mAb, unconjugated, detector	Human,Rat	ELISA
beta Amyloid	RB20120UC	beta Amyloid 1-38 Rabbit mAb, capture	Human,Mouse,Rat	ELISA
	RB20120UD	beta Amyloid 1-38 Rabbit mAb, unconjugated, detector	Human,Mouse,Rat	ELISA
	RB20121UC	beta Amyloid 1-42 Rabbit mAb, capture	Human,Mouse,Rat	ELISA
	RB20121UD	beta Amyloid 1-42 Rabbit mAb, unconjugated, detector	Human,Mouse,Rat	ELISA
HMGB1	MS20029UC	HMGB1 Mouse mAb, capture	Human,Mouse,Rat	ELISA
	MS20029HD	HMGB1 Mouse mAb, HRP conjugated, detector	Human,Mouse,Rat	ELISA