




配对抗体

Matched Antibody Pairs

- ① 抗体对经过严格的筛选配对和样本验证
- ② 高一致性、高灵敏度和高特异性确保获得最佳结果
- ③ 可用于夹心法 ELISA 和其他免疫检测
- ④ 包含常见细胞因子、肿瘤和其他疾病关键标志物等



正能生物



立足中国成都，服务全球用户

集抗体和蛋白研发、生产、销售于一体的生物高科技企业


正能官网 www.zen-bio.cn

成都正能生物（Zen-bioscience）成立于 2007 年，于 2013 年获得国家高新技术企业认证。公司拥有先进的蛋白研究方法，建成了集工具类蛋白及抗体研发生产为一体的全流程平台，囊括多个核心技术。目前已获得 8 项国家发明专利，34 项实用新型专利，并有多项发明专利已进入实审。产品覆盖了现代生命科学研究的各大前沿领域，包括癌症、免疫学、神经科学、心血管疾病、肿瘤干细胞、表观遗传学和内分泌代谢等。

正能生物致力于为全球生命科学研究机构及生物技术企业研制高品质的蛋白及抗体产品，从而促进生命科学和人类健康事业的发展！




ELISA检测原理示意图



配对抗体 (Matched Antibody Pair) 指可以同时结合在一个抗原分子上识别不同表位的两个抗体，一般应用于双抗夹心法 ELISA，流式液相多重检测 (Cytometric Bead Array, CBA)，也可适用于侧向层析 (LFA) 或其他相关的免疫检测。

正能提供的抗体对为小鼠单克隆抗体或重组兔单克隆抗体配对，均经过 ELISA 检测验证。抗体对产品线覆盖细胞因子、免疫球蛋白、肿瘤标志物和其他疾病相关蛋白等。我们还可提供匹配的蛋白质标准品。

部分产品的 ELISA 验证图



产品信息说明

货号前缀

MS：抗体来源小鼠
RB：抗体来源为兔
RT：抗体来源为大鼠

货号后缀

UC：未标记的捕获抗体
UD：未标记的检测抗体
HD：带 HRP 标记的检测抗体

产品名称后缀

capture：捕获抗体
detector：检测抗体

货号	产品名称	反应种属	应用
MS20010UC	AMBP Mouse mAb, capture	Human	ELISA
MS20010UD	AMBP Mouse mAb, unconjugated, detector	Human	ELISA
MS20011UC	CGB3 Mouse mAb, capture	Human	ELISA
MS20011UD	CGB3 Mouse mAb, unconjugated, detector	Human	ELISA
MS20012UC	CTLA4 Mouse mAb, capture	Human	ELISA
MS20012UD	CTLA4 Mouse mAb, unconjugated, detector	Human	ELISA
MS20013UC	Endostatin Mouse mAb, capture	Human	ELISA
MS20013UD	Endostatin Mouse mAb, unconjugated, detector	Human	ELISA
MS20014UC	Factor B Mouse mAb, capture	Human	ELISA
MS20014UD	Factor B Mouse mAb, unconjugated, detector	Human	ELISA
MS20015UC	FGF2 Mouse mAb, capture	Human	ELISA
MS20015UD	FGF2 Mouse mAb, unconjugated, detector	Human	ELISA
RB20015UC	FGF2 Rabbit mAb, capture	Human	ELISA
RB20015UD	FGF2 Rabbit mAb, unconjugated, detector	Human	ELISA
MS20016UC	G-CSF Mouse mAb, capture	Human	ELISA
MS20016UD	G-CSF Mouse mAb, unconjugated, detector	Human	ELISA
RB20016UC	G-CSF Rabbit mAb, capture	Human	ELISA
RB20016UD	G-CSF Rabbit mAb, unconjugated, detector	Human	ELISA
RT20017UC	GM-CSF Rat mAb, capture	Human	ELISA
RT20017UD	GM-CSF Rat mAb, unconjugated, detector	Human	ELISA
RB20017UC	GM-CSF Rabbit mAb, capture	Human	ELISA
RB20017UD	GM-CSF Rabbit mAb, unconjugated, detector	Human	ELISA
MS20018UC	IgA Mouse mAb, capture	Human	ELISA
MS20018UD	IgA Mouse mAb, unconjugated, detector	Human	ELISA
MS20019UC	IgE Mouse mAb, capture	Human	ELISA
MS20019UD	IgE Mouse mAb, unconjugated, detector	Human	ELISA
MS20020UC	IL-11 Mouse mAb, capture	Human	ELISA
MS20020UD	IL-11 Mouse mAb, unconjugated, detector	Human	ELISA
MS20021UC	IL-6 Mouse mAb, capture	Human	ELISA
MS20021UD	IL-6 Mouse mAb, unconjugated, detector	Human	ELISA
RB20021UC	IL-6 Rabbit mAb, capture	Human	ELISA
RB20021UD	IL-6 Rabbit mAb, unconjugated, detector	Human	ELISA
RT20022UC	Kit ligand Rat mAb, capture	Human	ELISA
RT20022UD	Kit ligand Rat mAb, unconjugated, detector	Human	ELISA
MS20023UC	Lactoferrin Mouse mAb, capture	Human	ELISA
MS20023UD	Lactoferrin Mouse mAb, unconjugated, detector	Human	ELISA
MS20024UC	IgM Mouse mAb, capture	Rat	ELISA
MS20024UD	IgM Mouse mAb, unconjugated, detector	Rat	ELISA
MS20025UC	Kappa light chain Mouse mAb, capture	Rat	ELISA
MS20025UD	Kappa light chain Mouse mAb, unconjugated, detector	Rat	ELISA
MS20026UC	IFN gamma Mouse mAb, capture	Human	ELISA
MS20026UD	IFN gamma Mouse mAb, unconjugated, detector	Human	ELISA
MS20027UC	IL-2 Mouse mAb, capture	Human	ELISA
MS20027UD	IL-2 Mouse mAb, unconjugated, detector	Human	ELISA
MS20028UC	TNF alpha Mouse mAb, capture	Human	ELISA
MS20028UD	TNF alpha Mouse mAb, unconjugated, detector	Human	ELISA
RB20028UC	TNF alpha Rabbit mAb, capture	Human	ELISA
RB20028UD	TNF alpha Rabbit mAb, unconjugated, detector	Human	ELISA
MS20030UC	CEA Mouse mAb, capture	Human	ELISA
MS20030UD	CEA Mouse mAb, unconjugated, detector	Human	ELISA
MS20031UC	ENO1 Mouse mAb, capture	Human	ELISA
MS20031UD	ENO1 Mouse mAb, unconjugated, detector	Human	ELISA
MS20032UC	PSA Mouse mAb, capture	Human	ELISA
MS20032UD	PSA Mouse mAb, unconjugated, detector	Human	ELISA
MS20033UC	MMP2 Mouse mAb, capture	Human	ELISA
MS20033UD	MMP2 Mouse mAb, unconjugated, detector	Human	ELISA
MS20036UC	S100B Mouse mAb, capture	Human	ELISA
MS20036UD	S100B Mouse mAb, unconjugated, detector	Human	ELISA
RT20034UC	Growth hormone 1 Rat mAb, capture	Human	ELISA
RT20034UD	Growth hormone 1 Rat mAb, unconjugated, detector	Human	ELISA
MS20035UC	EGF Mouse mAb, capture	Human	ELISA
MS20035UD	EGF Mouse mAb, unconjugated, detector	Human	ELISA
RT20037UC	Fibrinogen alpha Chain Rat mAb, capture	Human	ELISA
RT20037UD	Fibrinogen alpha Chain Rat mAb, unconjugated, detector	Human	ELISA
MS20038UC	alpha 1 Fetoprotein Mouse mAb, capture	Human	ELISA
MS20038UD	alpha 1 Fetoprotein Mouse Mouse mAb, unconjugated, detector	Human	ELISA
MS20039UC	IL-8 Mouse mAb, capture	Human	ELISA
MS20039UD	IL-8 Mouse mAb, unconjugated, detector	Human	ELISA
MS20040UC	IL-1 beta Mouse mAb, capture	Human	ELISA
MS20040UD	IL-1 beta Mouse mAb, unconjugated, detector	Human	ELISA
RB20040UC	IL-1 beta Rabbit mAb, capture	Human	ELISA
RB20040UD	IL-1 beta Rabbit mAb, unconjugated, detector	Human	ELISA
MS20041UC	Kappa Light Chain Mouse mAb, capture	Human	ELISA
MS20041UD	Kappa Light Chain Mouse mAb, unconjugated, detector	Human	ELISA
MS20042UC	IgG Mouse mAb, capture	Human	ELISA
MS20042UD	IgG Mouse mAb, unconjugated, detector	Human	ELISA

货号	产品名称	反应种属	应用
MS20043UC	IgM Mouse mAb, capture	Human	ELISA
MS20043UD	IgM Mouse mAb, unconjugated, detector	Human	ELISA
MS20044UC	MCP-1 Mouse mAb, capture	Human	ELISA
MS20044UD	MCP-1 Mouse mAb, unconjugated, detector	Human	ELISA
RB20044UC	MCP-1 Rabbit mAb, capture	Human	ELISA
RB20044UD	MCP-1 Rabbit mAb, unconjugated, detector	Human	ELISA
MS20046UC	Cystatin C Mouse mAb, capture	Human	ELISA
MS20046UD	Cystatin C Mouse mAb, unconjugated, detector	Human	ELISA
MS20047UC	FABP3 Mouse mAb, capture	Human	ELISA
MS20047UD	FABP3 Mouse mAb, unconjugated, detector	Human	ELISA
MS20048UC	EGF Mouse mAb, capture	Human	ELISA
MS20048UD	EGF Mouse mAb, unconjugated, detector	Human	ELISA
MS20049UC	Lipocalin 2 Mouse mAb, capture	Human,Rat	ELISA
MS20049UD	Lipocalin 2 Mouse mAb, unconjugated, detector	Human,Rat	ELISA
MS20050UC	RBP4 Mouse mAb, capture	Human	ELISA
MS20050UD	RBP4 Mouse mAb, unconjugated, detector	Human	ELISA
MS20051UC	Transferrin Mouse mAb, capture	Human	ELISA
MS20051UD	Transferrin Mouse mAb, unconjugated, detector	Human	ELISA
MS20052UC	FABP1 Mouse mAb, capture	Human	ELISA
MS20052UD	FABP1 Mouse mAb, unconjugated, detector	Human	ELISA
MS20053UC	Natriuretic Peptides B Mouse mAb, capture	Human	ELISA
MS20053UD	Natriuretic Peptides B Mouse mAb, unconjugated, detector	Human	ELISA
MS20054UC	Tau Mouse mAb, capture	Human	ELISA
MS20054UD	Phospho-Tau (Thr181) Mouse mAb, unconjugated, detector	Human	ELISA
MS20055UC	Tau Mouse mAb, capture	Human	ELISA
MS20055UD	Phospho-Tau (Thr217) Mouse mAb, unconjugated, detector	Human	ELISA
MS20056UC	Tau Mouse mAb, capture	Human	ELISA
MS20056UD	Phospho-Tau (Thr231) Mouse mAb, unconjugated, detector	Human	ELISA
RT20057UC	IgG1 Rat mAb, capture	Human	ELISA
MS20057UD	IgG Mouse mAb, unconjugated, detector	Human	ELISA
MS20058UC	IgG2 Mouse mAb, capture	Human	ELISA
MS20058UD	IgG2 Mouse mAb, unconjugated, detector	Human	ELISA
MS20059UC	IgG3 Mouse mAb, capture	Human	ELISA
MS20059UD	IgG Mouse mAb, unconjugated, detector	Human	ELISA
MS20060UC	IgG4 Mouse mAb, capture	Human	ELISA
MS20060UD	IgG4 Mouse mAb, unconjugated, detector	Human	ELISA
RB20061UC	IL-8 Rabbit mAb, capture	Human	ELISA
RB20061UD	IL-8 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20062UC	Cytokeratin 19 Rabbit mAb, capture	Human	ELISA
RB20062UD	Cytokeratin 19 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20063UC	CCL18 Rabbit mAb, capture	Human	ELISA
RB20063UD	CCL18 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20064UC	GRO alpha Rabbit mAb, capture	Human	ELISA
RB20064UD	GRO alpha Rabbit mAb, unconjugated, detector	Human	ELISA
RB20065UC	Thrombomodulin Rabbit mAb, capture	Human	ELISA
RB20065UD	Thrombomodulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20066UC	CXCL9 Rabbit mAb, capture	Human	ELISA
RB20066UD	CXCL9 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20068UC	IFN alpha1 Rabbit mAb, capture	Human	ELISA
RB20068UD	IFN alpha1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20069UC	IL-1 alpha Rabbit mAb, capture	Human	ELISA
RB20069UD	IL-1 alpha Rabbit mAb, unconjugated, detector	Human	ELISA
RB20070UC	IL-10 Rabbit mAb, capture	Human	ELISA
RB20070UD	IL-10 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20071UC	IL-12 p70/IL-12A Rabbit mAb, capture	Human	ELISA
RB20071UD	IL-12 p70/IL-12A Rabbit mAb, unconjugated, detector	Human	ELISA
RB20072UC	IL-12B Rabbit mAb, capture	Human	ELISA
RB20072UD	IL-12B Rabbit mAb, unconjugated, detector	Human	ELISA
RB20073UC	IL-13 Rabbit mAb, capture	Human	ELISA
RB20073UD	IL-13 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20074UC	IL-15 Rabbit mAb, capture	Human	ELISA
RB20074UD	IL-15 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20075UC	IL-16 Rabbit mAb, capture	Human	ELISA
RB20075UD	IL-16 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20076UC	IL-17A Rabbit mAb, capture	Human	ELISA
RB20076UD	IL-17A Rabbit mAb, unconjugated, detector	Human	ELISA
RB20077UC	IL-17F Rabbit mAb, capture	Human	ELISA
RB20077UD	IL-17F Rabbit mAb, unconjugated, detector	Human	ELISA
RB20078UC	IL-18 Rabbit mAb, capture	Human	ELISA
RB20078UD	IL-18 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20079UC	IL-1RA Rabbit mAb, capture	Human	ELISA
RB20079UD	IL-1RA Rabbit mAb, unconjugated, detector	Human	ELISA
RB20080UC	IL-22 Rabbit mAb, capture	Human	ELISA
RB20080UD	IL-22 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20081UC	IL-2RA Rabbit mAb, capture	Human	ELISA
RB20081UD	IL-2RA Rabbit mAb, unconjugated, detector	Human	ELISA

货号	产品名称	反应种属	应用
RB20082UC	IL-3 Rabbit mAb, capture	Human	ELISA
RB20082UD	IL-3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20083UC	IL-4 Rabbit mAb, capture	Human	ELISA
RB20083UD	IL-4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20084UC	IL-5 Rabbit mAb, capture	Human	ELISA
RB20084UD	IL-5 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20085UC	IL-7 Rabbit mAb, capture	Human	ELISA
RB20085UD	IL-7 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20087UC	MMP9 Rabbit mAb, capture	Human	ELISA
RB20087UD	MMP9 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20088UC	IL-9 Rabbit mAb, capture	Human	ELISA
RB20088UD	IL-9 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20089UC	IFN alpha2 Rabbit mAb, capture	Human	ELISA
RB20089UD	IFN alpha2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20090UC	LIF Rabbit mAb, capture	Human	ELISA
RB20090UD	LIF Rabbit mAb, unconjugated, detector	Human	ELISA
RB20091UC	M-CSF Rabbit mAb, capture	Human	ELISA
RB20091UD	M-CSF Rabbit mAb, unconjugated, detector	Human	ELISA
RB20092UC	MIF Rabbit mAb, capture	Human	ELISA
RB20092UD	MIF Rabbit mAb, unconjugated, detector	Human	ELISA
RB20093UC	MIP-1 alpha Rabbit mAb, capture	Human	ELISA
RB20093UD	MIP-1 alpha Rabbit mAb, unconjugated, detector	Human	ELISA
RB20094UC	RANTES Rabbit mAb, capture	Human	ELISA
RB20094UD	RANTES Rabbit mAb, unconjugated, detector	Human	ELISA
RB20095UC	TNF beta Rabbit mAb, capture	Human	ELISA
RB20095UD	TNF beta Rabbit mAb, unconjugated, detector	Human	ELISA
RB20096UC	IL-1 beta Rabbit mAb, capture	Mouse	ELISA
RB20096UD	IL-1 beta Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20097UC	IL-10 Rabbit mAb, capture	Mouse	ELISA
RB20097UD	IL-10 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20098UC	IL-12 p40 Rabbit mAb, capture	Mouse	ELISA
RB20098UD	IL-12 p40 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20099UC	IL-13 Rabbit mAb, capture	Mouse	ELISA
RB20099UD	IL-13 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20100UC	IL-17A Rabbit mAb, capture	Mouse	ELISA
RB20100UD	IL-17A Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20101UC	IL-2 Rabbit mAb, capture	Mouse	ELISA
RB20101UD	IL-2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20102UC	IL-4 Rabbit mAb, capture	Mouse	ELISA
RB20102UD	IL-4 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20103UC	IL-5 Rabbit mAb, capture	Mouse	ELISA
RB20103UD	IL-5 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20104UC	IL-6 Rabbit mAb, capture	Mouse	ELISA
RB20104UD	IL-6 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20105UC	IL-8 Rabbit mAb, capture	Mouse	ELISA
RB20105UD	IL-8 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20106UC	IFN gamma Rabbit mAb, capture	Mouse	ELISA
RB20106UD	IFN gamma Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20107UC	TNF alpha Rabbit mAb, capture	Mouse	ELISA
RB20107UD	TNF alpha Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20108UC	VEGFA Rabbit mAb, capture	Mouse	ELISA
RB20108UD	VEGFA Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20109UC	IL-1 beta Rabbit mAb, capture	Rat	ELISA
RB20109UD	IL-1 beta Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20110UC	IL-10 Rabbit mAb, capture	Rat	ELISA
RB20110UD	IL-10 Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20111UC	IL-6 Rabbit mAb, capture	Rat	ELISA
RB20111UD	IL-6 Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20112UC	TNF alpha Rabbit mAb, capture	Rat	ELISA
RB20112UD	TNF alpha Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20113UC	MIP-1 beta Rabbit mAb, capture	Human	ELISA
RB20113UD	MIP-1 beta Rabbit mAb, unconjugated, detector	Human	ELISA
RB20115UC	EpCAM Rabbit mAb, capture	Human	ELISA
RB20115UD	EpCAM Rabbit mAb, unconjugated, detector	Human	ELISA
RB20116UC	Alpha-synuclein Rabbit mAb, capture	Human	ELISA
RB20116UD	Alpha-synuclein Rabbit mAb, unconjugated, detector	Human	ELISA
RB20117UC	Alpha-synuclein Rabbit mAb, capture	Mouse	ELISA
RB20117UD	Alpha-synuclein Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20118UC	Angiopoietin 2/ANG2 Rabbit mAb, capture	Mouse	ELISA
RB20118UD	Angiopoietin 2/ANG2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20119UC	BACE1 Rabbit mAb, capture	Human	ELISA
RB20119UD	BACE1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20120UC	beta Amyloid 1-38 Rabbit mAb, capture	Human,Mouse,Rat	ELISA
RB20120UD	beta Amyloid 1-38 Rabbit mAb, unconjugated, detector	Human,Mouse,Rat	ELISA
RB20121UC	beta Amyloid 1-42 Rabbit mAb, capture	Human,Mouse,Rat	ELISA
RB20121UD	beta Amyloid 1-42 Rabbit mAb, unconjugated, detector	Human,Mouse,Rat	ELISA

货号	产品名称	反应种属	应用
RB20122UC	beta III Tubulin Rabbit mAb, capture	Human	ELISA
RB20122UD	beta III Tubulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20123UC	beta Tubulin Rabbit mAb, capture	Human	ELISA
RB20123UD	beta Tubulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20124UC	CD147 Rabbit mAb, capture	Human	ELISA
RB20124UD	CD147 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20125UC	CD147 Rabbit mAb, capture	Mouse	ELISA
RB20125UD	CD147 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20126UC	CD40 Rabbit mAb, capture	Human	ELISA
RB20126UD	CD40 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20127UC	CD40 Rabbit mAb, capture	Mouse	ELISA
RB20127UD	CD40 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20128UC	CD68 Rabbit mAb, capture	Human	ELISA
RB20128UD	CD68 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20129UC	COMP/Cartilage oligomeric matrix protein Rabbit mAb, capture	Human	ELISA
RB20129UD	COMP/Cartilage oligomeric matrix protein Rabbit mAb, unconjugated, detector	Human	ELISA
RB20130UC	Doublecortin Rabbit mAb, capture	Human	ELISA
RB20130UD	Doublecortin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20131UC	E Cadherin Rabbit mAb, capture	Mouse	ELISA
RB20131UD	E Cadherin Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20132UC	E Cadherin Rabbit mAb, capture	Rat	ELISA
RB20132UD	E Cadherin Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20133UC	E Cadherin Rabbit mAb, capture	Human	ELISA
RB20133UD	E Cadherin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20134UC	FOXA2 Rabbit mAb, capture	Human	ELISA
RB20134UD	FOXA2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20135UC	GAD65 Rabbit mAb, capture	Human	ELISA
RB20135UD	GAD65 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20136UC	GAD67 Rabbit mAb, capture	Human	ELISA
RB20136UD	GAD67 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20137UC	gamma Synuclein/SNCG Rabbit mAb, capture	Human	ELISA
RB20137UD	gamma Synuclein/SNCG Rabbit mAb, unconjugated, detector	Human	ELISA
RB20138UC	GAP43 Rabbit mAb, capture	Human	ELISA
RB20138UD	GAP43 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20139UC	GFAP Rabbit mAb, capture	Human	ELISA
RB20139UD	GFAP Rabbit mAb, unconjugated, detector	Human	ELISA
RB20140UC	GFAP Rabbit mAb, capture	Mouse	ELISA
RB20140UD	GFAP Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20141UC	MIF Rabbit mAb, capture	Mouse	ELISA
RB20141UD	MIF Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20142UC	MMP1 Rabbit mAb, capture	Human	ELISA
RB20142UD	MMP1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20143UC	MMP10 Rabbit mAb, capture	Mouse	ELISA
RB20143UD	MMP10 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20144UC	MMP12 Rabbit mAb, capture	Mouse	ELISA
RB20144UD	MMP12 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20145UC	MMP12 Rabbit mAb, capture	Human	ELISA
RB20145UD	MMP12 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20146UC	MMP13 Rabbit mAb, capture	Human	ELISA
RB20146UD	MMP13 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20147UC	MMP2 Rabbit mAb, capture	Mouse	ELISA
RB20147UD	MMP2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20148UC	MMP2 Rabbit mAb, capture	Human	ELISA
RB20148UD	MMP2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20150UC	MMP3 Rabbit mAb, capture	Mouse	ELISA
RB20150UD	MMP3 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20151UC	MMP3 Rabbit mAb, capture	Human	ELISA
RB20151UD	MMP3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20152UC	MMP3 Rabbit mAb, capture	Rat	ELISA
RB20152UD	MMP3 Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20153UC	MMP8 Rabbit mAb, capture	Human	ELISA
RB20153UD	MMP8 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20154UC	MMP8 Rabbit mAb, capture	Mouse	ELISA
RB20154UD	MMP8 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20155UC	MMP9 Rabbit mAb, capture	Human	ELISA
RB20155UD	MMP9 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20156UC	MMP9 Rabbit mAb, capture	Mouse	ELISA
RB20156UD	MMP9 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20157UC	Myelin oligodendrocyte glycoprotein Rabbit mAb, capture	Mouse	ELISA
RB20157UD	Myelin oligodendrocyte glycoprotein Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20158UC	N Cadherin Rabbit mAb, capture	Mouse	ELISA
RB20158UD	N Cadherin Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20159UC	Nestin Rabbit mAb, capture	Mouse	ELISA
RB20159UD	Nestin Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20160UC	Neurofilament Rabbit mAb, capture	Human	ELISA
RB20160UD	Neurofilament Rabbit mAb, unconjugated, detector	Human	ELISA

货号	产品名称	反应种属	应用
RB20161UC	Neurogranin Rabbit mAb, capture	Human	ELISA
RB20161UD	Neurogranin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20162UC	NG2 Rabbit mAb, capture	Mouse	ELISA
RB20162UD	NG2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20163UC	NG2 Rabbit mAb, capture	Human	ELISA
RB20163UD	NG2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20164UC	NMDAR1 Rabbit mAb, capture	Mouse	ELISA
RB20164UD	NMDAR1 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20165UC	NPTX2 Rabbit mAb, capture	Human	ELISA
RB20165UD	NPTX2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20166UC	NSE Rabbit mAb, capture	Human	ELISA
RB20166UD	NSE Rabbit mAb, unconjugated, detector	Human	ELISA
RB20167UC	NSE Rabbit mAb, capture	Mouse	ELISA
RB20167UD	NSE Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20168UC	P75 Neurogranin Rabbit mAb, capture	Human	ELISA
RB20168UD	P75 Neurogranin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20169UC	p75 NGF Receptor Rabbit mAb, capture	Mouse	ELISA
RB20169UD	p75 NGF Receptor Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20170UC	PDGFR alpha Rabbit mAb, capture	Human	ELISA
RB20170UD	PDGFR alpha Rabbit mAb, unconjugated, detector	Human	ELISA
RB20171UC	PDGFR beta Rabbit mAb, capture	Human	ELISA
RB20171UD	PDGFR beta Rabbit mAb, unconjugated, detector	Human	ELISA
RB20172UC	PDGFR beta Rabbit mAb, capture	Mouse	ELISA
RB20172UD	PDGFR beta Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20173UC	PGP9.5 Rabbit mAb, capture	Human	ELISA
RB20173UD	PGP9.5 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20174UC	PGP9.5 Rabbit mAb, capture	Mouse	ELISA
RB20174UD	PGP9.5 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20175UC	RAGE Rabbit mAb, capture	Rat	ELISA
RB20175UD	RAGE Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20176UC	RAGE Rabbit mAb, capture	Mouse	ELISA
RB20176UD	RAGE Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20177UC	S100B Rabbit mAb, capture	Human	ELISA
RB20177UD	S100B Rabbit mAb, unconjugated, detector	Human	ELISA
RB20178UC	S100A1 Rabbit mAb, capture	Human	ELISA
RB20178UD	S100A1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20179UC	S100A12/CGRP Rabbit mAb, capture	Human	ELISA
RB20179UD	S100A12/CGRP Rabbit mAb, unconjugated, detector	Human	ELISA
RB20180UC	S100A4 Rabbit mAb, capture	Human	ELISA
RB20180UD	S100A4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20181UC	S100A4 Rabbit mAb, capture	Mouse	ELISA
RB20181UD	S100A4 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20182UC	S100A9 Rabbit mAb, capture	Human	ELISA
RB20182UD	S100A9 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20183UC	SOX2 Rabbit mAb, capture	Human	ELISA
RB20183UD	SOX2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20184UC	SOX9 Rabbit mAb, capture	Human	ELISA
RB20184UD	SOX9 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20185UC	Stathmin 1 Rabbit mAb, capture	Human	ELISA
RB20185UD	Stathmin 1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20186UC	Tau Rabbit mAb, capture	Human	ELISA
RB20186UD	Tau Rabbit mAb, unconjugated, detector	Human	ELISA
RB20187UC	Tau Rabbit mAb, capture	Mouse	ELISA
RB20187UD	Tau Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20188UC	TDP43 Rabbit mAb, capture	Human	ELISA
RB20188UD	TDP43 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20189UC	TREM1 Rabbit mAb, capture	Human	ELISA
RB20189UD	TREM1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20190UC	TREM2 Rabbit mAb, capture	Mouse	ELISA
RB20190UD	TREM2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20191UC	TREM2 Rabbit mAb, capture	Human	ELISA
RB20191UD	TREM2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20192UC	VEGF Receptor 1 Rabbit mAb, capture	Human	ELISA
RB20192UD	VEGF Receptor 1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20193UC	VEGF Receptor 1 Rabbit mAb, capture	Rat	ELISA
RB20193UD	VEGF Receptor 1 Rabbit mAb, unconjugated, detector	Rat	ELISA
RB20194UC	VEGF Receptor 2 Rabbit mAb, capture	Human	ELISA
RB20194UD	VEGF Receptor 2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20195UC	VEGF Receptor 2 Rabbit mAb, capture	Mouse	ELISA
RB20195UD	VEGF Receptor 2 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20196UC	VEGF Receptor 3 Rabbit mAb, capture	Mouse	ELISA
RB20196UD	VEGF Receptor 3 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20197UC	VEGF Receptor 3 Rabbit mAb, capture	Human	ELISA
RB20197UD	VEGF Receptor 3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20198UC	VEGFB Rabbit mAb, capture	Mouse	ELISA
RB20198UD	VEGFB Rabbit mAb, unconjugated, detector	Mouse	ELISA

货号	产品名称	反应种属	应用
RB20199UC	Vimentin Rabbit mAb, capture	Human	ELISA
RB20199UD	Vimentin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20200UC	YKL-40/CHI3L1 Rabbit mAb, capture	Mouse	ELISA
RB20200UD	YKL-40/CHI3L1 Rabbit mAb, unconjugated, detector	Mouse	ELISA
RB20201UC	YKL-40/CHI3L1 Rabbit mAb, capture	Human	ELISA
RB20201UD	YKL-40/CHI3L1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20209UC	Insulin Rabbit mAb, capture	Human	ELISA
RB20209UD	Insulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20208UC	IL-20 Rabbit mAb, capture	Human	ELISA
RB20208UD	IL-20 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20218UC	Amyloid Precursor Protein Rabbit mAb,capture	Human	ELISA
RB20218UD	Amyloid Precursor Protein Rabbit mAb,unconjugated, detector	Human	ELISA
RB20219UC	Phospho-Tau (Thr181) Rabbit mAb,capture	Human	ELISA
RB20219UD	Phospho-Tau (Thr181) Rabbit mAb,unconjugated, detector	Human	ELISA
RB20220UC	Phospho-Tau (Thr231) Rabbit mAb,capture	Human	ELISA
RB20220UD	Phospho-Tau (Thr231) Rabbit mAb,unconjugated, detector	Human	ELISA
RB20221UC	Phospho-Tau (Ser202) Rabbit mAb,capture	Human	ELISA
RB20221UD	Phospho-Tau (Ser202) Rabbit mAb,unconjugated, detector	Human	ELISA
RB20222UC	IL-6 Rabbit mAb,capture	Human	ELISA
RB20222UD	IL-6 Rabbit mAb,unconjugated, unconjugated, detector	Human	ELISA
RB20223UC	ErbB2/HER2 Rabbit mAb,capture	Human	ELISA
RB20223UD	ErbB2/HER2 Rabbit mAb,unconjugated, detector	Human	ELISA
RT20063UC	IL-10 Rat mAb, capture	Human	ELISA
RT20063UD	IL-10 Rat mAb, detector	Human	ELISA
RB20010UC	Alpha-1-Microglobulin Rabbit mAb, capture	Human	ELISA
RB20010UD	Alpha-1-Microglobulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20027UC	IL-2 Rabbit mAb, capture	Human	ELISA
RB20027UD	IL-2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20038UC	alpha 1 Fetoprotein Rabbit mAb, capture	Human	ELISA
RB20038UD	alpha 1 Fetoprotein Rabbit mAb, unconjugated, detector	Human	ELISA
RB20049UC	Lipocalin-2/NGAL Rabbit mAb, capture	Human	ELISA
RB20049UD	Lipocalin-2/NGAL Rabbit mAb, unconjugated, detector	Human	ELISA
RB20050UC	RBP4 Rabbit mAb, capture	Human	ELISA
RB20050UD	RBP4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20034UC	Growth hormone Rabbit mAb, capture	Human	ELISA
RB20034UD	Growth hormone Rabbit mAb, unconjugated, detector	Human	ELISA
RB20224UC	SFN Rabbit mAb, capture	Human	ELISA
RB20224UD	SFN Rabbit mAb, unconjugated, detector	Human	ELISA
RB20225UC	Alkaline Phosphatase Rabbit mAb, capture	Human	ELISA
RB20225UD	Alkaline Phosphatase Rabbit mAb, unconjugated, detector	Human	ELISA
RB20226UC	beta 2 Microglobulin Rabbit mAb, capture	Human	ELISA
RB20226UD	beta 2 Microglobulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20227UC	C Reactive Protein Rabbit mAb, capture	Human	ELISA
RB20227UD	C Reactive Protein Rabbit mAb, unconjugated, detector	Human	ELISA
RB20228UC	Calcitonin Rabbit mAb, capture	Human	ELISA
RB20228UD	Calcitonin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20229UC	CD14 Rabbit mAb, capture	Human	ELISA
RB20229UD	CD14 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20230UC	CD4 Rabbit mAb, capture	Human	ELISA
RB20230UD	CD4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20231UC	CD55 Rabbit mAb, capture	Human	ELISA
RB20231UD	CD55 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20232UC	CEACAM1 Rabbit mAb, capture	Human	ELISA
RB20232UD	CEACAM1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20233UC	CEACAM6 Rabbit mAb, capture	Human	ELISA
RB20233UD	CEACAM6 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20234UC	CTHRC1 Rabbit mAb, capture	Human	ELISA
RB20234UD	CTHRC1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20235UC	Cytokeratin 18 Rabbit mAb, capture	Human	ELISA
RB20235UD	Cytokeratin 18 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20236UC	beta Amyloid 1-40 Rabbit mAb, capture	Human	ELISA
RB20236UD	beta Amyloid 1-40 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20237UC	DPP4/CD26 Rabbit mAb, capture	Human	ELISA
RB20237UD	DPP4/CD26 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20238UC	Elafin/ESI Rabbit mAb, capture	Human	ELISA
RB20238UD	Elafin/ESI Rabbit mAb, unconjugated, detector	Human	ELISA
RB20239UC	Eotaxin 2 Rabbit mAb, capture	Human	ELISA
RB20239UD	Eotaxin 2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20240UC	Eph receptor A2 Rabbit mAb, capture	Human	ELISA
RB20240UD	Eph receptor A2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20241UC	EPO Rabbit mAb, capture	Human	ELISA
RB20241UD	EPO Rabbit mAb, unconjugated, detector	Human	ELISA
RB20242UC	FGF21 Rabbit mAb, capture	Human	ELISA
RB20242UD	FGF21 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20243UC	Fibrinogen Rabbit mAb, capture	Human	ELISA
RB20243UD	Fibrinogen Rabbit mAb, unconjugated, detector	Human	ELISA

货号	产品名称	反应种属	应用
RB20244UC	Fibulin 1 Rabbit mAb, capture	Human	ELISA
RB20244UD	Fibulin 1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20245UC	Glypican 1/ GPC1 Rabbit mAb, capture	Human	ELISA
RB20245UD	Glypican 1/ GPC1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20246UC	HE4 Rabbit mAb, capture	Human	ELISA
RB20246UD	HE4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20247UC	HGF Rabbit mAb, capture	Human	ELISA
RB20247UD	HGF Rabbit mAb, unconjugated, detector	Human	ELISA
RB20248UC	HLA-DR Rabbit mAb, capture	Human	ELISA
RB20248UD	HLA-DR Rabbit mAb, unconjugated, detector	Human	ELISA
RB20250UC	ICAM2 Rabbit mAb, capture	Human	ELISA
RB20250UD	ICAM2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20251UC	ICAM3 Rabbit mAb, capture	Human	ELISA
RB20251UD	ICAM3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20252UC	IDH1 Rabbit mAb, capture	Human	ELISA
RB20252UD	IDH1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20253UC	IGFBP2 Rabbit mAb, capture	Human	ELISA
RB20253UD	IGFBP2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20254UC	IGFBP7 Rabbit mAb, capture	Human	ELISA
RB20254UD	IGFBP7 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20255UC	IL-17B Rabbit mAb, capture	Human	ELISA
RB20255UD	IL-17B Rabbit mAb, unconjugated, detector	Human	ELISA
RB20256UC	IL-17C Rabbit mAb, capture	Human	ELISA
RB20256UD	IL-17C Rabbit mAb, unconjugated, detector	Human	ELISA
RB20257UC	IL-17D Rabbit mAb, capture	Human	ELISA
RB20257UD	IL-17D Rabbit mAb, unconjugated, detector	Human	ELISA
RB20259UC	IL-25 Rabbit mAb, capture	Human	ELISA
RB20259UD	IL-25 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20260UC	IL-33 Rabbit mAb, capture	Human	ELISA
RB20260UD	IL-33 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20261UC	TIM 1 Rabbit mAb, capture	Human	ELISA
RB20261UD	TIM 1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20262UC	L1CAM Rabbit mAb, capture	Human	ELISA
RB20262UD	L1CAM Rabbit mAb, unconjugated, detector	Human	ELISA
RB20263UC	LAMC2 Rabbit mAb, capture	Human	ELISA
RB20263UD	LAMC2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20264UC	LECT2 Rabbit mAb, capture	Human	ELISA
RB20264UD	LECT2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20265UC	Ki67 Rabbit mAb, capture	Human	ELISA
RB20265UD	Ki67 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20266UC	Liver FABP Rabbit mAb, capture	Human	ELISA
RB20266UD	Liver FABP Rabbit mAb, unconjugated, detector	Human	ELISA
RB20267UC	LOXL2 Rabbit mAb, capture	Human	ELISA
RB20267UD	LOXL2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20268UC	Midkine Rabbit mAb, capture	Human	ELISA
RB20268UD	Midkine Rabbit mAb, unconjugated, detector	Human	ELISA
RB20269UC	MUC1 Rabbit mAb, capture	Human	ELISA
RB20269UD	MUC1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20270UC	MUC16 Rabbit mAb, capture	Human	ELISA
RB20270UD	MUC16 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20271UC	Neutrophil Elastase Rabbit mAb, capture	Human	ELISA
RB20271UD	Neutrophil Elastase Rabbit mAb, unconjugated, detector	Human	ELISA
RB20272UC	OLFM4 Rabbit mAb, capture	Human	ELISA
RB20272UD	OLFM4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20273UC	Osteocalcin Rabbit mAb, capture	Human	ELISA
RB20273UD	Osteocalcin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20274UC	Osteopontin Rabbit mAb, capture	Human	ELISA
RB20274UD	Osteopontin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20275UC	Parathyroid Hormone Rabbit mAb, capture	Human	ELISA
RB20275UD	Parathyroid Hormone Rabbit mAb, unconjugated, detector	Human	ELISA
RB20276UC	Pentraxin 3/PTX3 Rabbit mAb, capture	Human	ELISA
RB20276UD	Pentraxin 3/PTX3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20277UC	Pepsinogen II/PGC Rabbit mAb, capture	Human	ELISA
RB20277UD	Pepsinogen II/PGC Rabbit mAb, unconjugated, detector	Human	ELISA
RB20278UC	PICP Rabbit mAb, capture	Human	ELISA
RB20278UD	PICP Rabbit mAb, unconjugated, detector	Human	ELISA
RB20279UC	Plexin B1 Rabbit mAb, capture	Human	ELISA
RB20279UD	Plexin B1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20280UC	PR3 Rabbit mAb, capture	Human	ELISA
RB20280UD	PR3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20281UC	Pro-GRP Rabbit mAb, capture	Human	ELISA
RB20281UD	Pro-GRP Rabbit mAb, unconjugated, detector	Human	ELISA
RB20282UC	Procalcitonin Rabbit mAb, capture	Human	ELISA
RB20282UD	Procalcitonin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20283UC	Procollagen Type 1 Rabbit mAb, capture	Human	ELISA
RB20283UD	Procollagen Type 1 Rabbit mAb, unconjugated, detector	Human	ELISA

货号	产品名称	反应种属	应用
RB20284UC	Progranulin Rabbit mAb, capture	Human	ELISA
RB20284UD	Progranulin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20285UC	Prostate Specific Antigen Rabbit mAb, capture	Human	ELISA
RB20285UD	Prostate Specific Antigen Rabbit mAb, unconjugated, detector	Human	ELISA
RB20286UC	Ribonuclease 3/ECP Rabbit mAb, capture	Human	ELISA
RB20286UD	Ribonuclease 3/ECP Rabbit mAb, unconjugated, detector	Human	ELISA
RB20288UC	SerpinB3/SCCA Rabbit mAb, capture	Human	ELISA
RB20288UD	SerpinB3/SCCA Rabbit mAb, unconjugated, detector	Human	ELISA
RB20289UC	SPARC Rabbit mAb, capture	Human	ELISA
RB20289UD	SPARC Rabbit mAb, unconjugated, detector	Human	ELISA
RB20290UC	Survivin Rabbit mAb, capture	Human	ELISA
RB20290UD	Survivin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20291UC	Tartrate Resistant Acid Phosphatase Rabbit mAb, capture	Human	ELISA
RB20291UD	Tartrate Resistant Acid Phosphatase Rabbit mAb, unconjugated, detector	Human	ELISA
RB20292UC	Tenascin C Rabbit mAb, capture	Human	ELISA
RB20292UD	Tenascin C Rabbit mAb, unconjugated, detector	Human	ELISA
RB20293UC	TFPI Rabbit mAb, capture	Human	ELISA
RB20293UD	TFPI Rabbit mAb, unconjugated, detector	Human	ELISA
RB20294UC	TGF beta 3 Rabbit mAb, capture	Human	ELISA
RB20294UD	TGF beta 3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20295UC	TGFBI Rabbit mAb, capture	Human	ELISA
RB20295UD	TGFBI Rabbit mAb, unconjugated, detector	Human	ELISA
RB20296UC	Thrombospondin 2 Rabbit mAb, capture	Human	ELISA
RB20296UD	Thrombospondin 2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20297UC	TIMP2 Rabbit mAb, capture	Human	ELISA
RB20297UD	TIMP2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20298UC	uPA Receptor/U-PAR Rabbit mAb, capture	Human	ELISA
RB20298UD	uPA Receptor/U-PAR Rabbit mAb, unconjugated, detector	Human	ELISA
RB20299UC	VCAM1 Rabbit mAb, capture	Human	ELISA
RB20299UD	VCAM1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20300UC	Visfatin Rabbit mAb, capture	Human	ELISA
RB20300UD	Visfatin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20301UC	Vitamin D Binding protein Rabbit mAb, capture	Human	ELISA
RB20301UD	Vitamin D Binding protein Rabbit mAb, unconjugated, detector	Human	ELISA
RB20302UC	IGFBP4 Rabbit mAb, capture	Human	ELISA
RB20302UD	IGFBP4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20303UC	DKK4 Rabbit mAb, capture	Human	ELISA
RB20303UD	DKK4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20202UC	IL-23 Rabbit mAb, capture	Human	ELISA
RB20202UD	IL-23 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20311UC	GAPDH Rabbit mAb, capture	Human	ELISA
RB20311UD	GAPDH Rabbit mAb, unconjugated, detector	Human	ELISA
RB20312UC	Thymidine Kinase 1/TK1 Rabbit mAb, capture	Human	ELISA
RB20312UD	Thymidine Kinase 1/TK1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20313UC	ESM1 Rabbit mAb, capture	Human	ELISA
RB20313UD	ESM1 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20314UC	Granzyme B Rabbit mAb, capture	Human	ELISA
RB20314UD	Granzyme B Rabbit mAb, unconjugated, detector	Human	ELISA
RB20315UC	HAVCR2 Rabbit mAb, capture	Human	ELISA
RB20315UD	HAVCR2 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20316UC	GNLY Rabbit mAb, capture	Human	ELISA
RB20316UD	GNLY Rabbit mAb, unconjugated, detector	Human	ELISA
RB20317UC	CCL20 Rabbit mAb, capture	Human	ELISA
RB20317UD	CCL20 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20318UC	LAG-3 Rabbit mAb, capture	Human	ELISA
RB20318UD	LAG-3 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20319UC	Granzyme A Rabbit mAb, capture	Human	ELISA
RB20319UD	Granzyme A Rabbit mAb, unconjugated, detector	Human	ELISA
RB20320UC	CIQA Rabbit mAb, capture	Human	ELISA
RB20320UD	CIQA Rabbit mAb, unconjugated, detector	Human	ELISA
RB20321UC	CTLA4 Rabbit mAb, capture	Human	ELISA
RB20321UD	CTLA4 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20322UC	CXCL13 Rabbit mAb, capture	Human	ELISA
RB20322UD	CXCL13 Rabbit mAb, unconjugated, detector	Human	ELISA
RB20323UC	Glycodelin Rabbit mAb, capture	Human	ELISA
RB20323UD	Glycodelin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20324UC	CEA Rabbit mAb, capture	Human	ELISA
RB20324UD	CEA Rabbit mAb, unconjugated, detector	Human	ELISA
RB20325UC	Albumin Rabbit mAb, capture	Human	ELISA
RB20325UD	Albumin Rabbit mAb, unconjugated, detector	Human	ELISA
RB20001UC	Proinsulin Rabbit mAb, capture	Human	ELISA
RB20001UD	Proinsulin Rabbit mAb, unconjugated, detector	Human	ELISA



PlatinumAb
精品抗体



RecRabs®
重组免单抗



Histo-plus®
病理级IHC抗体



Antibody pairs
抗体对



KO/KD validated
KO/KD验证抗体



mIHC kit
多重免疫荧光

成都正能生物技术有限公司



扫码关注我们

☎ 4008 863973

🌐 www.zen-bio.cn


地 址：四川省·成都市·高新区·天府生命科技园

邮 箱：support@zen-bioscience.com (技术支持)

sales@zen-bioscience.com (订购咨询)



正能各区域代理商联系方式



扫码查看代理商

①登录“正能官网”→ 顶部菜单栏【关于我们】→【联系我们】

②关注“正能生物”公众号 → 底部菜单栏【抗体订购】→【订购渠道】

③拨打“咨询热线”，转1

